


theradome



We understand how difficult hair loss can be

Theradome® proudly restores the hair and confidence of people across the globe. *We can help you too.*



Which stage of hair loss are you in?

Your stage of hair loss will help determine which regimen will help you to achieve *your hair goal*


Hair loss Factors




Many factors contribute to hair loss

- Stress
- Poor nutrition
- Hormones
- Genetic factors
- Environmental factors
- Inflammation

FEMALE PATTERN BALDNESS



MALE PATTERN BALDNESS







theradome

Your hair health goals are unique to you. *What are yours?*


Reverse the process
of thinning -
miniaturization



Minimize hair loss -
reverse the process
of shedding



Experience new
hair growth



What are your treatment options?



LIGHT
THERAPY



PRP
INJECTIONS



TOPICAL
TREATMENTS



INSTANT
VOLUMIZING
SOLUTIONS



PRESCRIPTION
MEDICATION




SURGICAL
INTERVENTION

Low Level Light Therapy (LLLT) is the only Clinical Proven, Natural, Drug-Free, Pain Free, Non-Invasive hair growth treatment option with no side effects and cost effective in the long term.



LASER PHOTOTHERAPY

Red laser light penetrates deep in the scalp and reaches hair follicle roots with the proper wavelength at 680nm




LASER THERAPY boosts your results

- ✓ Hair shedding decreases
- ✓ Hair thickness increases
- ✓ Scalp health improves
- ✓ New hair grows

Stimulating your follicles with the right wavelength of red visible light will do these three things:

1. Slow or stop hair loss
2. Thicken your existing hair follicles
3. Help you regrow new hair






See results with the
Theradome® PRO LH80
in as soon as 4 weeks!

- ✓ Flexible + convenient
- ✓ Simple to use
- ✓ Hands-free
- ✓ Clinically tested
- ✓ At home use
- ✓ Light weight

Theradome® PRO LH80 hair restoration device

Technologically Advanced Laser Hair Growth Treatment for safe and effective home use. Clinical trials have shown that hair loss can be reversed with the Theradome PRO LH80 Laser Helmet. Result can be seen in as soon as 4 weeks although this may vary from person to person.






See results with the
 Theradome® EVO LH40
 in as soon as 4 weeks!

- ✓ Flexible + convenient
- ✓ Simple to use
- ✓ Hands-free
- ✓ Clinically tested
- ✓ At home use
- ✓ Light weight

Theradome® EVO LH40 hair restoration device

Technologically Advanced Laser Hair Growth Treatment for safe and effective home use. Clinical trials have shown that hair loss can be reversed with the Theradome EVO LH40 Laser Helmet. Result can be seen in as soon as 4 weeks although this may vary from person to person.






Why is Theradome® different?

quality
— OVER —
~~quantity~~

Theradome® delivers the most powerful OTC device treatment available.


Theradome's proprietary VL680 Laser is the deepest penetrating FDA-cleared laser available to those suffering from hair loss. The fact is, when more joules are transferred to the scalp, hair follicles receive greater growth stimulation.





Why is Theradome[®] different?

The less divergent a beam of light energy, the more focused it is and the greater penetration and absorption it will have. The VL680 is 2.4x less divergent and therefore more focused, which allows for better penetration.



What does this mean? *Theradome[®] delivers stronger results faster*

IT REALLY WORKS! *Real people, real results*

98%

of Theradome® users experience positive results



- Reverse miniaturization
- Minimize hair loss
- Grow new hair



"Best Investment I have ever made. Amazing product, THANK YOU THERADOME."

Glenn Hill, Los Angeles, CA



"I am so thankful for this product. It has made a huge difference in my part size. Thank you Theradome for giving me my hair back!"

Teresa Brown, Woodside, NY



"The front of my scalp has started to regrow hair covering up an area that was becoming very noticeable."

J.D., Denver, CO



"More hair growth everywhere. Love it!!"

J.S., Franklin, TN



"This is easier to use than the hand held version and my hair looks thicker. Wish I had this 5 years ago!"

Aaron Goodwill, Christiansted, VI

GROW NEW
HAIR WITH
THERADOME®
2 X 20 MINUTES
PER WEEK



START YOUR
JOURNEY

theradome

Day
1



Day
180




I'm very grateful for Theradome. The back of my head is really filling in with new hair. I have more confidence now, it really works!

James Di Paolo, AZ



THERADOME®
DELIVERS
AWARD
WINNING
TECHNOLOGY,
SAFE AND
EFFECTIVE FOR



START YOUR
JOURNEY

theradome

Day
1



Day
180




Theradome has made my hair thicker and fuller in only 6 months. I would definitely recommend it for hair regrowth.

Joann Whitefield, TX




GROW NEW
HAIR WITH
THERADOME®
2 X 20 MINUTES
PER WEEK



START YOUR
JOURNEY

theradome

Day
1



Day
120



Theradome has already in a short amount of time made new hair grow that wasn't there before.

Mary Aguilera, VA





The team at Theradome® are
passionate about hair growth

We value the people that partner with us on
their journey towards thicker, healthier hair

Join us today!



theradome
www.theradome.in

+91-8800712312, +91-11-40109681 Info@theradome.in

AS SEEN IN

ELLE

Women's Health

EONLINE

SKIN,INC
Advanced Treatment Science

abc

SHAPE

CBS

LIFETIME


FOX
BUSINESS

The New York Times

BAZAR


hba

theradome



Doctor recommended Laser Hair
Growth Treatment for Safe and
Effective Home Use

 Info@theradome.in

 +91-8800712312, +91-11-40109681

 www.theradome.in

 162, Ansal Chamber -II, Bhikaji
Cama Place New Delhi- 110066

